



PATENT  
Attorney Docket No.: 27373/39055B

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Applicant: Wang *et al.*

Application No. 09/748,710

Filed: December 22, 2000

For: METHOD FOR  
GENERATION OF LONGER  
CDNA FRAGMENTS FROM  
SAGE TAGS FOR GENE  
IDENTIFICATION

Group Art Unit: 1637

Examiner: Joyce Tung

) I hereby certify that this correspondence is  
) being deposited with the U.S. Postal Service  
) with sufficient postage as First Class Mail, in  
) an envelope addressed to: Commissioner for  
) Patents, P.O. Box 1450, Alexandria, VA  
) 22313-1450, on the date shown below.

) Date: August 26, 2003

) William K. Muhl

**DECLARATION OF SAN MING WANG, PH.D., UNDER 37 C.F.R. §1.131**

I, San Ming Wang, Ph.D., state that,

1. I am an inventor of the subject matter claimed in the above-identified application, which claims priority to U.S. Provisional Application Nos.: 60/174,391 and 60/173,617, filed January 3, 2000, and December 29, 1999, respectively. I am making this declaration to provide evidence that the subject matter claimed in the above-identified application was completed in the United States at a date at least prior to July 19, 1999, which is the date of acceptance of van den Berg, et al., Nucleic Acids Research, 1999, Vol. 27(17), pages i-iii, cited by the Examiner in an Office Action mailed March 26, 2003.

2. To establish the date of completion of the invention, copies of pages from my laboratory notebook are attached as Exhibit A. All dates appearing on these laboratory notebook pages have been redacted.

3. The pages from my laboratory notebook establish the following facts:

(a) Information from SAGE tags alone do not provide sufficient information for gene identification;

(b) The solution to the problems associated with SAGE tags was to extend the SAGE tag sequence to the 3' end of the cDNA;


(c) The longer 3' cDNA could be generated using a SAGE tag as a primer in an extension and/or amplification reaction(s);

(d) One amplification reaction useful in solving the SAGE tag problem is PCR, using a SAGE tag as a sense primer and a set of universal primers consisting of oligo (dT) with A, G, CA, CG, CC anchors as antisense primers, and using the cDNA sample used in SAGE tag collection as a substrate;


4. These laboratory notebook pages confirm that the subject matter of the pending claims in the application was completed at least prior to July 19, 1999.

5. This declaration is submitted in response to a non-final Office Action.

6. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true and, further, that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.



San Ming Wang, M.D.



Date



**EXHIBIT A**



select candidate genes from different classes based on the TAG sequences from Lin Zhang, John Hopkins.

• Confirm the matched genes from tag seq. by Blastn parameter set: expect: 1000, cutoff: 60.

sequences aligned show many aligned genes in many cases to define the most likely ones, it needs to stress:


a. 5' add CATG

e.g. <sup>5'</sup> CATG XXXXXXXXXXXX <sub>tag</sub>

b. the matched sequences should be located in the last N/A site of the sequence <sup>3'</sup>

c. confirm the right gene by <sup>generally</sup> using large templates based on the tag sequence. This is not available in SAGE technique. I designed a system to stress this issue. If solved, it will largely improve the reality of ~~seq~~ gene identification in SAGE.

REDACTED



- based on Tag seq. design 5' seq sense primer, plus 6 bases tail to maintain the stringency


~ anchored dT  
dy. It is known well  
+ Tag pol  
precisely distinguish  
base anchored primer  
CR reactions

- use anchored primer as 3' antisense primer to select 3' subpopulation for PCR. IMPORTANT! NEW IDEA.
- PCR amplify cDNA
- the cDNA should be pre-digested by Nla III. The possible outcome would be
- specific temp amplified because of the 5' sense primer / 3' anchored primer gives exponential amplification, others will only be ~~double~~ increased, i.e. 3. e. - potential problem ~;
- Amplification efficiency lower due to 3' extension consume large portions of substrates
- specificity ?


REDACTED

model system to confirm the possibility.

genes  $\beta$ -Actin 52/AS ~~S2/Act~~ ~~S2/GAT~~  
HSC 70 mixture of AT, A/G/C  
HSP 75  
HSP 78





antisense 0.5 (AS) 0.5 (mix arch act) ~ ~ ~ ~ ~



95 20" 95 20" } x 25  
50 20" } x 5 = 60 20" } x 25  
72 20" 72 20"

conclusion. the mixture of a-chained primers doesn't provide  
and ...

The possible reason for failure to ~~to~~ amplify may be that the mixture of  $\alpha$ -tubulin disturbed the efficient amplification; Try to use separated  $\alpha$ -tubulin as 3' primer, use ACT1 as model



$\beta$ -actin mRNA end (HSAC07)

TAAGTG CACAC CTTA TTTTT  
 A TTTTT (GACT1)

Conclusion (1)  $\alpha$ -tubulin as 3' C.K.

improve:

- increase 3' anchor primer amount
- decrease first 5 cycles in PCR to favor 3' primer bind to temp, then shift to 60°C for annealing
- increase cycle number to 40 ✓ Actin


	8x
Myc/2	40
buffer	40
QNTIP	40
Taq	8
CDNA	8
Sense	8

DATA  
HT<sub>in</sub>  
ACT<sub>in</sub>

ef1-  
ribosomal  
(3) x  
S24

X037K7  
M69013  
U26726

AA5142  
A219316  
S28  
AA67692  
AA61797  
AJ223473  
D25786



antiseise 1 1 5 10 15 15 10 15 1 5 10 15  
50ng 250ng 500ng 750ng

420 14 14 10 5 - 14 10 5 - 14 10 5 -  
25.2L per

94°C 20"  
42°C 20"  
72°C 20" } x5 → 60°C 20" } x30  
72°C 20"

S/As ATn GTn CTn  
50ng 250 500 750






Conclusion:

- ① increase anchor primer amount did increase the amplification eff.
  - ② CTn still generated diffused
- in principle, this strategy should work! Can be use




to turn it into ~~single~~ <sup>double</sup> lines, and print it as  
the other template appears at least, and  
perhaps better, as follows:



\* Also, the cDNA should be generated with anchor set, not oligo dT which creates background in PCR.

It seems that one base anchored dT primers may not generate enough specificity in particular temp, due to the complicated cDNA pool composition, or different members in the same family like Hsp70s.

To increase the specificity, it may be useful to use 2 base anchored dT, set as 3' antisense primer. The complete set of sense / antisense will be



- To identify each ~~temp~~ tag gene, the total reaction will be 11 PCR reaction with Tag seq as Sense (+ 6 base 5' to 3' direction).
- If negative, reversed Tag sequence will be used as the Sense primer in case the original tag was generated in 3'-5' direction.
- In this way, the specific 3' fragment will be generated, and match to database.

next. try repeat exp. with anchor primer generated cDNA → { one base anchor 5' 5' }